

# Techniques for BSC Decontamination

Henry S. Luftman, PhD

DRS Laboratories, Inc.

1-888-377-1533

[hluftman@drslaboratories.com](mailto:hluftman@drslaboratories.com)

# Outline

- General steps in performing a BSC Decontamination
- Chemical-Specific Issues
  - Methodologies
  - Advantages / Disadvantages
- Comparison
- DRS Laboratories' Position

# Choice of Decontaminants

- Formaldehyde Gas
- Hydrogen Peroxide Vapor
- Chlorine Dioxide Gas
- Others
  - Methyl Bromide
  - Ethylene Oxide
  - Ozone

# Typical Applications

- Clean-up of a contamination event
- Maintenance or other need to access contaminated plenums
- Moving of BSC
- HEPA filter replacement
- End of work program

# Requirements for a Successful SD

- Decontamination – typically looking for a log 4-6 reduction of test bacterial spores
- Choice of decontaminant
- Penetration to all surfaces
- Penetration through HEPA filter and into “dead legs”
- Temperature and humidity control
- Containment of fumigant

# Requirements for a Successful SD (Cont.)

- Disposal of decontaminator
  - Vent, neutralize, scrub
- Validation of decontamination
  - Biological indicators
- Material compatibility
- Safety

# General Preparation Options

- Seal BSC when decontamination required
- Construct BSC with decontamination facilities included
- Permanently modify BSC for specific decon type
- Ensure gas-tight damper if ducted to building
- Insert recirculation (optional)

# General Preparation – Biological Indicators

- Given 7 –day test time, typically assume decontamination method already validated
- If using indicators, often *B. atrophaeus* or *G. stearothermophilus*
- Use appropriate substrate (not cellulose for HP)
- Upstream and/or downstream of HEPA filters
- Log-Reduction enumeration vs. Go / No Go
- Controls



# General Preparation - Final

- Establish and measure proper humidity & temperature
- Final seal
- Pressure check (neutral to adjacent area)
- Establish safety perimeter
- Meet OSHA requirements

# General Procedure

- Fumigant generation to a steady state concentration
- Environmental monitoring for leakage
  - E.g., Draeger pumps, infrared analyzers
- Appropriate personnel protective equipment (PPE)
  - Full face respirator, gloves, lab coat
- Neutralization or scrubbing → Ventilation

# General Procedure (Cont.)

- Validation of BSC Decontamination
  - Biological Indicators (opt.)
  - Monitoring of relative humidity, space temperature, and/or decontaminant concentration during process

# Formaldehyde Gas (CH<sub>2</sub>O)

- Typically via depolymerization of Paraformaldehyde (PF)
- NSF standard 0.3 gm/ft<sup>3</sup> → ~8000 ppm
- Mechanism: methylation of DNA
- Requires relative humidity > 60%
- Target contact time > 6 hr
- Use *Bacillus atrophaeus* as BI

# Formaldehyde Gas (Cont.)

- Neutralization with ammonia gas ( $\text{NH}_3$ )
  - ~ 1 hr contact time
- Vent and environmental monitoring
- Clean “fall-out”
  - Mixture of methenamine and PF
  - Can limit PF with humidity control

# Formaldehyde Gas - Advantages

- “True” gas
- Relatively inexpensive
- General material compatibility
- Industry accepted

# Formaldehyde Gas - Issues

- “Fall-out” residue
  - Added clean-up time
- Carcinogen
- Potential odor residual
- Polymerization on cold surfaces

# Hydrogen Peroxide Vapor ( $\text{H}_2\text{O}_2$ )

- Typically delivered by flash vaporization of aqueous peroxide mixture
  - The mixture is generally close to or above saturation in air
- Two major vendors of generators with significant differences
  - Mechanism : Oxidation
- Required contact time less than formaldehyde
- Use *Geobacillus stearothermophilus* as BI



# HP Vapor – STERIS (VHP)

- Avoids condensation on surfaces to minimize corrosion and optimize distribution
- Typically two portals into BSC for VHP inlet and return
- Design cabinet with appropriate circulation paths
- Dehumidify to  $< 30\%$  RH

# HP Vapor – STERIS (Cont.)

- Typical 1-2mg/liter, 750-1500 ppm (D~1-2 min)
  - Data that D value is lower than for liquid HP
- Target 70-85% RH during decontamination
- Continually introduce HP, decomposing HP in return
- Cycle Phases
  - Dehumidification / Conditioning / Decontamination / Aeration

# HP Vapor – BIOQUELL (Clarus)

- Seeks “micro-condensation”
  - BQ believes D ~2 min required liquid presence
- Swiveling source to inject high-speed droplets to all surfaces
  - Condensate “bounces”
- Monitor for onset of condensation

# HP Vapor Advantages

- Safe by-products (water and oxygen)
- No residue
- Industry accepted
- Automated
- Relatively short cycle time if properly engineered

# HP Vapor Issues

- Instability of HP toward decomposition
- Decomposition may block access of decontaminant
- Condensation may cause control issues
- Cellulose materials absorb or decompose
  - May effect decontamination or aeration
- Some material issues – nylon, cellulose, copper, lead, iron oxide, epoxy
  - Condensation may effect painted surfaces
- Capital equipment cost

# HP Operating Conditions

- (Based on limited information)
- Typical duration from conditioning through decontamination:
  - 1-2 hours
- Typical aeration:
  - 2-4 hours

# Chlorine Dioxide Gas (ClO<sub>2</sub>)

- Mechanism: Selective oxidation (no chloridation)
- Generated on site via reaction
  - $\text{Cl}_2(\text{g}) + 2\text{NaClO}_2 \longrightarrow 2\text{ClO}_2(\text{g}) + 2\text{NaCl}$
- Visible green gas
- Humidification required, 65-90% RH
- D-value 0.1-0.8 min for 10-30 mg/L
  - (3,500-10,000 ppm)

# Chlorine Dioxide Gas (Cont.)

- Scrubbing
  - Wet, with alkaline solutions
  - Dry, via absorption (e.g., charcoal)
- Direct venting option (use by paper industry)
- Monitor concentration and relative humidity
- Use *Geobacillus stearothermophilus* as BI (?)
- The Halide Group, ClorDiSys Solutions, Sabre Technologies (different reaction)



# Chlorine Dioxide Gas – Advantages

- Safe by-products (oxygen and salt)
- No residue
- Not flammable / explosive
- “True gas – no condensation issues
- Reputation for use in Anthrax decontamination

# Chlorine Dioxide Gas - Issues

- Less well-known or characterized
- Mild corrosion/discoloring to cold steel, copper, brass
  - Particularly in the presence of water
- Potentially corrosive if Chlorine gas ( $\text{Cl}_2$ ) is present
  - Care to avoid  $\text{Cl}_2$  in synthesized CD
  - Care to avoid  $\text{Cl}_2$  creation by UV exposure
- Current low PEL limit (0.1 ppm)

# CD Operating Conditions

- Typical duration conditioning through decontamination:
  - 1-1.5 hr
- Typical duration of aeration / scrubbing:
  - 0.25-0.5 hr

# Comparison

Issue	Formaldehyde Gas	Hydrogen Peroxide Vapor	Chlorine Dioxide Gas
Sporocidal effectiveness	+	+	+
Effective through HEPA filters	+	?	+
Non Carcinogenic	-	+	+
Toxicity (TWA PEL)	0.75 ppm	1.0 ppm	0.1 ppm
Non-explosive (at used concentrations)	-	-	+
Humidity requirement (RH)	60-90%	30% (Steris) or ambient (Bioquell)	65-90%
No residue	-	+	+
Non-corrosive	+	+ (dry) / ? (cond.)	+ / - (with chlorine)
Method of removal	Neutralizer	Catalytic breakdown	Scrubbing
Limited development effort	+	-	+
Limited cost	+	-	-

# NSF/ANSI 49 – 2002

## Annex G

### **Recommended microbiological Decontamination procedure**

*“Prior to decontamination with an alternative [note added: other than depolymerized paraformaldehyde] method (such as VHP), cycle parameters and validation of those parameters must be developed for each model and size of BSC.”*

# Baker Company in CleanRooms (Mar. 2001)

*“Existing cabinets can be modified in the field to accept hydrogen peroxide vapor generators... However, these alternatives are less than optimal because:*

- 1. Such cabinets are not designed for fast distribution of  $H_2O_2$  and poor vapor distribution will require long cycle times – considerably longer than for formaldehyde cycle. ...”*


# DRS Laboratories' Position

- Safety issue in maintaining BSC following Clarus decontamination
- Need validation data for decontamination covering
  - Variation of safety cabinet
  - Variation of HEPA filter thickness
  - Demonstration of decontamination on up and down stream sides of HEPA filter
  - Demonstration of decontamination if cabinet has a failed motor

# Chlorine Dioxide Alternative

- DRS Laboratories has experience with two delivery systems:
  - Direct  $\text{ClO}_2$  insertion into BSC
    - Continuous concentration monitoring
    - Wet “scrubbing” of gas after cycle
  - Tablet generation of  $\text{ClO}_2$  in water
    - “Bubbling” air to draw out gas
    - Dry “scrubbing” after cycle
    - No additional cost relative to PF





BI Test Locations

# Current CD Data

- Four runs with four BI's each (*G. Stearothermophilus*)
- 125 min. from humidification to opening
- Exhaust plenum wo/ recirculation:
  - 3.5-5.5 log reduction
- Exhaust plenum w/ recirculation:
  - 6.0-6.2 log reduction
- Of others,
  - 11/12 > 6.2 log reduction, 1/12 ~ 5.6 log reduction

# Chlorine Dioxide Validation

- Less of an issue, as  $\text{ClO}_2$  is “true” gas
- Have begun process at DRS Laboratories on multiple cabinets
- Propose to monitor with BI's on other cabinets on short term after DRS Laboratories internal validation complete
- DRS Laboratories' technicians would not require waiting for BI results during this phase